

SAMPLE SUBMISSION TO GENOMICS UNIT

1. Fill the Sequencing or Microarrays service request form and the sample quality control request form and e-mail it to us (genomics@crg.eu)

You will find the requests forms on Genomics Unit web page (“Services & Equipment” tab):
<http://www.crg.eu/en/programmes-groups/genomics-unit>

2. We prefer 1.5 ml tube for the samples and should have:

- **Name** of the sample (please, use the same sample name as in the service request form). Max 4 characters, capital letters. Write sample name on tube-cap (do NOT write any information on tube walls: we will stick our internal number there).
- Minimal volume **10 ul** (we need to have enough volume for samples QC and the sequencing or microarray service)

3. Time schedule for receiving samples at Genomics Unit if you come in person:

Monday 11h-13h

Wednesday 11h-13h

Friday 11h-13h

If this days/time are not suitable for you, please let us know so we can agree another day / time for you.

4. Alternatively, if you send samples via courier, there is not a specific time schedule, but it is important to inform us in advance.

In this case, please direct them to:

Genomics Unit
Centre de Regulació Genòmica (CRG)
PRBB Building
C/ Dr. Aiguader 88, 5th floor, lab 520
08003 Barcelona

Please include inside the package the service request forms in paper and also e-mail the request form (electronic format) to us (genomics@crg.eu).

It is important that you let us know that you have send some samples, so we can be aware that we should receive a package.

Please, contact **Genomics Unit staff** (genomics@crg.eu) if you have any question or need further information, thanks!